


EasyOne Extractor



Avantages

- Ease of use:** Fully automatic R&D machine. The user sets pressure, temperature, and extraction time, then starts the process.
- Precision and control:** Touchscreen interface for easy control and monitoring of temperature and pressure, up to **350 bar** and **95°C**. This allows quick testing of pressure/temperature combinations to compare selectivity and optimize extraction parameters.
- Sample preparation:** Enables static supercritical CO₂ extraction before chromatographic analysis such as GC-MS, GC-O, HPLC, or LC-MS. Useful for building aromatic profiles or comparing chemical signatures.
- Application versatility:** Suitable for proof-of-concept studies and screening of P/T extraction conditions. EasyOne can also be used for maceration, impregnation, thermal cycles, and solubility studies under supercritical conditions.


Applications

Solubility studies
Extraction
Maceration
Impregnation
Thermal cycling



Contact the TENNAD team to design a solution:
contact@tennad.com

EasyOne Extractor



Compact R&D System for Supercritical CO₂:

EasyOne is a compact, fully autonomous R&D extractor for supercritical CO₂ experiments. It offers a simple way to explore CO₂ extraction in teaching, analytical and research laboratories.

Precise Control for Fast Screening:

EasyOne automatically controls pressure, temperature and extraction time. Users can quickly screen P/T conditions, compare selectivity and identify suitable parameters for a given biomass.

Easy Operation and Laboratory-Scale Design:

A touchscreen HMI simplifies operation and parameter monitoring. EasyOne is available with 125 mL or 250 mL baskets, with an automatic cleaning cycle for easier routine use.

Wide Range of R&D Applications:

EasyOne supports R&D screening, solubility studies, static CO₂ extraction, static CO₂ treatment, maceration, CO₂ assisted impregnation, thermal cycling and feasibility testing. A co-solvent may also be used depending on the application.

Sample Preparation for Chromatographic Analysis:

EasyOne can prepare extracts before GC-MS, GC-FID, GC-O, HPLC or LC-MS analysis. It helps generate clean, concentrated and reproducible extracts for aroma profiling and chemical fingerprinting.

A Simple Tool to Compare and Validate CO₂ Processes:

EasyOne helps users move from first trials to validated extraction conditions through fast P/T screening, controlled operation and reproducible results within one compact system.

Automatic R&D extractor ideal for static CO₂ extraction, fast P/T screening, and sample preparation in laboratories and education.



Contact the TENNAD team to design a solution:
contact@tennad.com